

physio sensing



PHYSIOSENSING RANGE TECHNOLOGY
Balance | Feet Pressure Map | Posturography
Virtual Reality | Visual Biofeedback

powered by

sensingfu+ure

Contents


Abouts us	4
PhysioSensing Range Technology	6
Pressure Plate	8
Force Plate	10
Podo Software	12
Balance Software	14
Virtual Reality Libra VR Clinic	16
Plate Comparison	18
Software Comparison	19
Define your solution	21
Ask us	22

We design, develop and implement medical devices in a unique synergy between Engineering and Health.


Specialists in technology for force and pressure platforms with wide applications in physical and vestibular rehabilitation.


Software with powerful balance&pressure analysis tools, oriented to the health professional and the patient.

Over 10 years of experience and certified by the ISO 13485 standard for medical devices and the ISO 9001 standard for the quality management system.

 www.physiosensing.net
www.sensingfuture.pt

 info@sensingfuture.pt

 +351 239 404 234

 Instituto Pedro Nunes - Bloco E
Rua Pedro Nunes
3030-199 Coimbra, Portugal

sensingfu+ure
a greater step



version 1.0 revision date: 09.2020

Our numbers

15

Developed projects and technologies

31

Countries with our technology

10

Years of experience in the development of technological medical devices



PhysioSensing Range Technology


1. Balance

Balance control consists of controlling the body center of mass over its limits of stability. Clinical balance assessment can help assess fall risk and/or determine the underlying reasons for balance disorders. The benefits of using force plates in balance assessment comes from their ability to measure center of pressure (CoP).

2. Comfortable stance with the eyes closed




2. Feet Pressure Map



Foot Pressure Mapping is a method of measuring pressures on the surface of the foot during standing or walking. Static and dynamic baropodometric analysis are performed on a baropodometric platform, through which the pressure exerted by the feet from standstill and during walking is being measured.

3. Posturography

Computerized posturography systems utilize force platforms to measure the sway of a patient by determining the movements of the instantaneous Center of Pressure (CoP). The CoP data collected can be visualized through a statokinesi-gram and stabilogram.



4. Virtual Reality

In the recent decade, Virtual Reality (VR) has become generally accepted as a therapeutic tool for neurological patients. VR involves real-time simulation and interactions between sensory, motor and cognitive channels. VR can be set up to be strongly immersive, in that the environment appears real and three-dimensional. VR provides an ideal environment to study the balancing strategies.



5. Visual Biofeedback

Visual Biofeedback (VB) is a rehabilitation method that can be used during static balance training, offering the patient visual information on the position of the center of gravity within the range of stability as the patient stands on a plate. VB stimulates motivation, proprioceptive information to the patient, simulates the body movements and provides valuable information for the health professional.



Pressure Plate

Light and portable device.
The most accurate and cost
effective way for
balance&pressure
assessment and
training.

PhysioSensing allows you
to evaluate your clinical
practice and make it
objective and quantified
in a clinical report.



Technology

- 1 cm² sensor
- 1600 sensors
- Baropodometric plate

Balance evaluation
and training



Technical Specifications



CE Medical Device Class I according to directive 93/42/EEC

Type	Portable
Size	61 x 58 cm
Thickness	1 cm
Weight	4 kg
Active Surface	40 x 40 cm
Sensors number	1600
Sensor size	1 x 1 cm
Sensor type	Resistive
Sensor life time	more than 1 000 000 actuations
Maximum pressure (each sensor)	100 N/cm ²
Temperature range	from 0°C to 60°C
Connection/power	USB
Frequency	100 Hz ~100 acquisitions/seconds

PC Requirements

CPU Processor	i3, Quad-Core, with 2GHz
RAM	4 GB
USB Ports	2.0 or 3.0
Operating System	Windows 7, 8, 10 (32 or 64-bit)

Compatible Software

Physiosensing
Balance Software

Physiosensing
Podo Software

Virtual Reality
Libra VR Clinic

Accessories

Carry Bag



Foam



40 x 40 x 8 cm
18 Kg/m³ density

Force Plate

Normalized stabilometric platform for postural rehabilitation.

Fundamental technology to clinical practice of balance disorders.

Extreme precision of the CoP and high sensitivity level



Technology

- 3 load cells
- Stabilometric plate
- In accordance with French Association of Posturology.

Balance evaluation and training



Technical Specifications



CE Medical Device Class I according to directive 93/42/EEC

Size (Length x Width x Height)	53 x 46 x 3,5 cm
Thickness	1,2 cm
Weight	7,8 kg
Material	Aluminium AU4G
Maximal load	128 kg
Non linearity	< 0,2 %
Resolution	900 points/Kg
Hysteresis	< 0,2 %
Sampling rate	Adjustable from 5 Hz to 40 Hz
Analogic / Digital conversion	16 bits
Platform computer interface	USB
Power supply	USB cable

PC Requirements

CPU Processor	i3 with 2GHz
RAM	2 GB
USB Ports	2.0 or 3.0
Operating System	Windows 7, 8, 10 (32 or 64-bit)

Compatible Software

- Physiosensing Balance Software
- Virtual Reality Libra VR Clinic

Accessories

Foam



40 x 40 x 8 cm
18 Kg/m3 density

Dynamic Balancing Plate



Used on top of the platform to perform dynamic ML or AP Balance.

Podo Software


Static and dynamic baropodometric analysis

Static & Postural Analysis

Analyze your plantar pressure map distribution in a single (static) image or during a defined time.




- Static mapping with center of pressure
- Max, Average pressure calculations
- Weight distribution
- Measuring possibilities (length, angle, pressure, area)
- Exam comparasion
- Statokinesigram and Stabilogram




Dynamic Analysis


Analyze your plantar pressure distribution on the platform during a one or four gait cycle



Dynamic Analysis Report




Frame view



Dynamic Analysis view

- Frame view
- CPEI
- CoP line
- Force, Area, Pressure, velocity Graphs
- PEAK risk
- Gait Cycle Phases



Exam Comparasion

Compatible Plates





Pressure Plate

General Features

- Clinical Reports with normative Data
- Report
- Patient Management
- Database Export (Excel, PDF)

Balance Software

Powerful Balance
assessment & training tool



mCTSIB- Modified Clinical Test of Sensory Interaction on Balance & Romberg Test

This protocol allows the static balance measurement in four sensorial conditions: stable surface and open eyes; stable surface and closed eyes; surface instable and eyes open and surface instable and eyes closed.

Body Sway

Use Body Sway to create a personalized posturography. Define initial conditions and obtain CoP variation, ML and AP variables over time. It also includes more than 30 parameters derived from a posturographic examination including Fourier analyzes.

Limits of Stability LOS

Perhaps the most used protocol for balance assesment. This protocol quantifies the directional control and the maximum distance that the patient can reach with its center of pressure in 8 different directions.

Fall Risk

Use Fall Risk to measure the static balance in four conditions: comfortable stance with eyes open and eyes closed, narrow stance with eyes open and closed. After performing, the value of the sway velocity index appears and provides a fall risk prediction.

Rhythmic Weight Shift

Use the Rhythmic Weight Shift protocol to evaluate the transfer capacity of the center of pressure rhythmically in the sagittal and anteroposterior plane, at three different velocities.

Unilateral Stance

Use this protocol to measure the balance in four conditions: left foot lifted up with eyes open, left foot lifted up with eyes closed, right foot lifted up with eyes open and right foot lifted up with eyes closed.

Balance Error Scoring System

The BESS protocol allows the measurement of postural stability with eyes closed in three different positions on two types of surface (firm and unstable): two feet together, unipodal and tandem position.



Assessment

Balance improvement training

In the Sagittal and Anteroposterior exercises, the patient must reach the balance position in the sagittal plane or in the anteroposterior plane.



In the Spiderweb, Tunnel, Dots pattern and Bars pattern exercises the goal is to keep a balance position despite the visual stimulus.



In the exercises Route Bars, Route Spiral, Route Square and Route Maze the goal is to go through the routes within its boundaries with the center of pressure, following the red dots.



In the exercises Square, Circle, Eight and Spiral, the patient must reach all the dots disposed in the form indicated in the name.



In the Follow the Point and Moving Route exercises, the goal is to reach the moving red point and follow it within its tolerance margin.



Compatible Plates



Pressure Plate



Force Plate

General Features

- Clinical Reports with normative Data
- Report
- Patient Management
- Database Export (Excel, PDF)

Virtual Reality Libra VR Clinic

A high-end virtual reality based system, designed for vestibular, balance and oculomotor disorders. Therapy becomes controlled, high customized and trackable.

Immersive stimuli

Exciting and engaging virtual worlds train the patient brain

Enhanced therapy

Through digital technology


Monitored tracking

Allows for assessing and adjusting exercises

Evidence-based and approved

Scientifically backed to
meet therapists' needs

Clinic dashboard
provides customization
and reporting



Virtual Reality
to deliver immersive
virtual reality stimulation
and exercises including:



- › Smooth Pursuit
- › Saccadic
- › Optokinetic Nystagmus
- › VOR
- › VOR suppression
- › Supermarket effect and Visual Parallax
- › Vergence
- › Fixation
- › Seek and Find games

Optional
Platform for posturography and
balance games

Compatible Plates



Pressure Plate



Force Plate

Required hardware

Virtual Reality Glasses

Computer







CPU Processor	i5, Quad-Core, with 3.30 GHz
RAM	8 GB
USB Ports	3.0
Operating System	Windows 10 (64-bit)
Display	DisplayPort 1.2 or Mini DisplayPort 1.2
Graphics Card	NVIDIA GTX 1050Ti or better VR Ready graphics card

Plate Comparison




	Pressure Plate	Force Plate
Portability	● ● ● ● ●	● ● ● ● ●
Dimensions Length x Width x Height	61 x 58 mm x 1 cm	53 x 46 x 3,5 cm
Weight	4 Kg	7,8 Kg
Technology	1600 resistive sensors	3 load cells
Type	Baropodometric Plate	Stabilometric Plate
CoP accuracy	● ● ● ● ●	● ● ● ● ●
Feet pressure map	Yes	No
If you valorize more	Feet Pressure Map	CoP Path, velocity, Area; AP&ML components
Software compatibility	<div><div>🖥️ Podo Software</div><div>🖥️ Balance Software</div><div>🖥️ Virtual Reality Libra VR Clinic</div></div>	<div><div>🖥️ Balance Software</div><div>🖥️ Virtual Reality Libra VR Clinic</div></div>

Software Comparison

	Podo Software	Balance Software
 Balance Assessment Protocols (mCTSIB, LOS, RWS, WBS, US, FR, BESS, Sit to Stand, RT, Body Sway-Posturography)		✓
 Balance Training Exercises (Sagital & Anteroposterior, Visual Stimulus, Static Figures, Static Paths, Random Static, LOS Training, Moving target, Follow me)		✓
 Static Analysis * This feature only works with Pressure Plate	✓	✓ *
 Postural Analysis	✓	
 Dynamic Analysis	✓	
 Balance games		✓

Define your solution



Devices can be used together or separately.
Choose the best setting.

I	Pressure Plate + Podo Software		PODO
II	Pressure Plate + Balance Software		REHAB
III	Pressure Plate + Balance Software + Podo Software		REHAB
IV	Pressure Plate + Balance Software + Virtual Reality Libra VR Clinic		VERTIGO
V	Force Plate + Balance Software		VERTIGO
VI	Force Plate + Balance Software + Virtual Reality Libra VR Clinic		VERTIGO
VII	Virtual Reality Libra VR Clinic		VERTIGO

+ Additional equipment

Mobile Solution A
Suitable for laptop computer. Adjustable height.



Mobile Solution B
Suitable for laptop or desktop computer. Possibility to add a larger monitor. Adjustable height.



Static Solution A
Suitable for wall mounting. You can have a laptop or desktop computer and even adapt a larger monitor.





— Need help
choosing
the right
solution?

**ask
us**

quotation
demonstration
redirect to
our distributor

info@sensingfuture.pt

www.physiosensing.net

we ship worldwide



sensingfu+ure

a greater step

www.physiosensing.net

